

Product Name:	VWR-HC-FFP-20 76579-146
Description:	<p>VWR Catalog Number: 76579-146 VWR Standard Series Natural Refrigerant Flammable Storage Freezer 20 Cu. Ft. Capacity One solid swing door, right hinged Microprocessor temperature controller Digital temperature display with sample simulated probe Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: -15°C – -25°C Six fixed shelves Keyed door lock Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation Natural, Hydrocarbon Refrigerant (R290) Two year parts and labor warranty Utilizes ATEX and non-sparking components for safety Exterior dimensions: 30 3/4" W x 29" D x 74 5/8" H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost 115V, 60 Hz, 6 Amps, 1/2 HP ETL listed Unit meets NFPA guidelines 45 and 70 Shipping weight: 240 lbs.</p>