



Material	23510-540
Material description	Chloroform
Grade	ACS
CAS Number	67-66-3
Molecular formula	CHCl ₃
Molecular mass	119,38

Characteristics

Specifications

Assay (additives included)	Min. 99,8 %
Acetone and aldehyde	Passes test
Acid and chloride	Passes test
Free chlorine	Passes test
Presence of preservative	Passes test
Substances darkened by sulphuric acid	Passes test
Color (Pt-Co)	Max. 10 APHA
Residue after evaporation	Max. 0,001 %
Pb (Lead)	Max. 0,05 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada