



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Date of issue: 11/17/2014

Version 1.0

## SECTION 1. Identification

### Product identifier

Product number	PX1565
Product name	Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS
CAS-No.	7778-77-0

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Reagent for analysis
-----------------	----------------------

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   General Inquiries: +1-978-715-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)
---------	--

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

## SECTION 2. Hazards identification

### GHS-Labeling

Not a dangerous substance according to GHS.

### Other hazards

None known.

## SECTION 3. Composition/information on ingredients

Formula	$\text{KH}_2\text{PO}_4$	$\text{H}_2\text{KO}_4\text{P}$ (Hill)
Molar mass	136.08 g/mol	

Remarks	No hazardous ingredients according to the OSHA Hazard Communication Standard 29 CFR 1910.1200.
---------	---

## SECTION 4. First aid measures

### Description of first-aid measures

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number  
Product name

PX1565  
Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS

Version 1. 0

### *Inhalation*

After inhalation: fresh air.

### *Skin contact*

After skin contact: wash off with plenty of water. Remove contaminated clothing.

### *Eye contact*

After eye contact: rinse out with plenty of water.

### *Ingestion*

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

Never give anything by mouth to an unconscious person.

### **Most important symptoms and effects, both acute and delayed**

irritant effects, gastric pain, Diarrhea, Nausea, Vomiting, Discomfort

### **Indication of any immediate medical attention and special treatment needed**

No information available.

---

## **SECTION 5. Fire-fighting measures**

### **Extinguishing media**

#### *Suitable extinguishing media*

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### **Special hazards arising from the substance or mixture**

Not combustible.

Ambient fire may liberate hazardous vapors.

Fire may cause evolution of:

Oxides of phosphorus

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

In the event of fire, wear self-contained breathing apparatus.

#### *Further information*

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Avoid inhalation of dusts. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### **Environmental precautions**

Do not empty into drains.

### **Methods and materials for containment and cleaning up**

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number  
Product name

PX1565  
Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS

Version 1.0

Cover drains. Collect, bind, and pump off spills.  
Observe possible material restrictions (see sections 7 and 10).  
Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

---

### SECTION 7. Handling and storage

#### Precautions for safe handling

Observe label precautions.

#### Conditions for safe storage, including any incompatibilities

Tightly closed. Dry.

Store at room temperature.

---

### SECTION 8. Exposure controls/personal protection

#### Exposure limit(s)

Contains no substances with occupational exposure limit values.

#### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

#### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

#### *Hygiene measures*

Change contaminated clothing. Wash hands after working with substance.

#### *Eye/face protection*

Safety glasses

#### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

#### *Respiratory protection*

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

### SECTION 9. Physical and chemical properties

Physical state	solid
Color	colorless
Odor	odorless
Odor Threshold	Not applicable

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number PX1565 Version 1.0  
Product name Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS

---

pH	4.2 - 4.6 at 20 g/l 68 °F ( 20 °C)
Melting point	ca. 253 °C (decomposition)
Boiling point/boiling range	Not applicable, (decomposition)
Flash point	does not flash
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	Not applicable
Upper explosion limit	Not applicable
Vapor pressure	No information available.
Relative vapor density	No information available.
Density	2.34 g/cm <sup>3</sup> at 68 °F ( 20 °C)
Relative density	No information available.
Water solubility	222 g/l at 68 °F ( 20 °C)
Partition coefficient: n-octanol/water	No information available.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Oxidizing properties	none
Ignition temperature	not combustible
Bulk density	ca. 1,200 kg/m <sup>3</sup>

---

### SECTION 10. Stability and reactivity

#### Reactivity

See below

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number  
Product name

PX1565  
Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS

Version 1.0

### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

### Possibility of hazardous reactions

Violent reactions possible with:

Strong oxidizing agents, Bases, acids

### Conditions to avoid

Strong heating (decomposition).

### Incompatible materials

no information available

### Hazardous decomposition products

in the event of fire: See section 5.

---

## SECTION 11. Toxicological information

### Information on toxicological effects

#### *Likely route of exposure*

Eye contact, Skin contact, Ingestion

#### *Acute oral toxicity*

LDLO Rat: 4,640 mg/kg (RTECS)

Symptoms: Nausea, Vomiting, gastric pain, Diarrhea, Discomfort

#### *Acute dermal toxicity*

LD50 Rabbit: > 4,640 mg/kg  
(RTECS)

#### *Eye irritation*

slight irritation

#### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

#### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

#### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

### Carcinogenicity

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number  
Product name

PX1565  
Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS

Version 1.0

ACGIH by NTP.  
No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

### Further information

Hazardous properties cannot be excluded but are unlikely when the product is handled appropriately.  
Handle in accordance with good industrial hygiene and safety practice.

---

## SECTION 12. Ecological information

### Ecotoxicity

#### *Toxicity to fish*

LC0 Leuciscus idus (Golden orfe): ca. 900 mg/l; 48 h (IUCLID)

### Persistence and degradability

No information available.

### Bioaccumulative potential

No information available.

### Mobility in soil

No information available.

#### *Additional ecological information*

Depending on the concentration, phosphates may contribute to the eutrophication of water supplies.

Discharge into the environment must be avoided.

---

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

## SECTION 14. Transport information

### Land transport (DOT)

Not classified as dangerous in the meaning of transport regulations.

### Air transport (IATA)

Not classified as dangerous in the meaning of transport regulations.

### Sea transport (IMDG)

Not classified as dangerous in the meaning of transport regulations.

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1565

Version 1. 0

Product name

Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS

---

### SECTION 15. Regulatory information

#### United States of America

##### SARA 313

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

##### SARA 302

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

##### Clean Water Act

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

##### DEA List I

Not listed

##### DEA List II

Not listed

#### US State Regulations

##### Massachusetts Right To Know

Remarks

No components are subject to the Massachusetts Right to Know Act.

##### California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

##### Notification status

TSCA: All components of the product are listed in the TSCA-inventory.

DSL: All components of this product are on the Canadian DSL.

---

### SECTION 16. Other information

#### Training advice

Provide adequate information, instruction and training for operators.

#### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Date of issue: 11/17/2014

---

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1565

Version 1. 0

Product name

Potassium Phosphate, Monobasic [Potassium Dihydrogen Phosphate] GR ACS

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*