



Item Number PO158
Item Potassium Bicarbonate, Granular, USP
CAS Number 298-14-6
Molecular Formula KHCO₃
Molecular Weight 100.12

Assay (Dried Basis)..... 99.5-101.5%

MAXIMUM LIMITS

Loss on Drying..... 0.3 %
Normal Carbonate..... 2.5%
Heavy Metals (as Pb)..... 0.001%
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843