



Office 337.237.7123

Toll Free 800.396.7123

Fax 337.237.8712

Material Safety Data Sheet Sodium Thiosulfate Solutions 0.025N to 1N

Section 1 - Chemical Product and Company Identification

MSDS Name: Sodium Thiosulfate Solutions 0.025N to 1N

Synonyms: None

Company Identification:

C & G Containers, Inc.
152 Easy Street
Lafayette, LA 70506

For information, call: 337-237-7123

Emergency Number: 201-796-7100

For CHEMTREC assistance, call: 800-424-9300

For International CHEMTREC assistance, call: 703-527-3887

Section 2 - Composition, Information on Ingredients

CAS#	Chemical Name	Percent	EINECS/ELINCS
7732-18-5	Water	76-99.4	231-791-2
10102-17-7	Sodium Thiosulfate, pentahydrate	0.6-22.9	unlisted

Section 3 - Hazards Identification

EMERGENCY OVERVIEW

Appearance: colorless liquid.

Caution! May cause eye, skin, and respiratory tract irritation. This is expected to be a low hazard for usual industrial handling.

Target Organs: None.

Potential Health Effects

Eye: May cause mild eye irritation.

Skin: May cause skin irritation.