



Sodium Hydroxide, 50% Solution AR

Product No. 7705
Specifications current as of: Jun 8
2011

TEST	SPECIFICATION
Assay (NaOH)	49.5 - 52 %
Identification	Passes Test
Ammonium Hydroxide Precipitate	max. 0.010 %
Chloride (Cl)	max. 0.0025 %
Heavy Metals (as Ag)	max. 0.0015 %
Iron (Fe)	max. 0.0005 %
Nitrogen Compounds (as N)	max. 0.0005 %
Potassium (K)	max. 0.025 %
Phosphate (PO ₄)	max. 0.0005 %
Sulfate (SO ₄)	max. 0.0025 %
Sodium Carbonate (Na ₂ CO ₃)	max. 0.1 %
Product Information (not specifications):	
Appearance (clear, colorless liquid)	
For Laboratory, Research or Manufacturing Use	
Protect from air to avoid absorption of carbon dioxide.	