

Material Safety Data Sheet L-Ascorbic Acid

Section 1 - Chemical Product and Company Identification

MSDS Name: L-Ascorbic Acid

Synonyms: Ascorbic acid; 3-Keto-L-glucofuranolactone; Vitamin C; L-(+)-Ascorbic acid

Company Identification:

C & G Container, Inc.
152 easy Street
Lafayette, LA 70506

For information, call: (337) 237-7123

Emergency Number: 201-796-7100

For CHEMTREC assistance, call: 800-424-9300

For International CHEMTREC assistance, call: 703-527-3887

Section 2 - Composition, Information on Ingredients

CAS#	Chemical Name	Percent	EINECS/ELINCS
50-81-7	L-Ascorbic acid	>99	200-066-2

Hazard Symbols: None listed.

Risk Phrases: None listed.

Section 3 - Hazards Identification

EMERGENCY OVERVIEW

Appearance: white to yellow crystals.

Caution! Strong reducing agent. Fire and explosion risk in contact with oxidizing agents. May cause eye, skin, and respiratory tract irritation. Light sensitive. Air sensitive.