



## Streptavidin conjugated to Dylight® 650

Product Number D10-2212  
Lot# DY650001-13-003  
Amount 1.0 mg  
Concentration 1.0 mg/mL  
Store at 2-8°C

### Form/Shipping & Storage

Supplied as a liquid. Upon receipt, store at 2-8°C in the dark, do not freeze.

### Handling

We recommend that the investigator determine the appropriate working concentration for their specific application.

Avoid exposure to heat and light.


### Buffer

The product is in 100 mM sodium phosphate (pH 7.4), 100 mM sodium chloride, 2mM sodium azide as a preservative.

For research use only, not for diagnostic or therapeutic use.

### Spectral Characteristics

Visible absorption maxima 655  
Emission maximum 672  
F/P Ratio 2.84



**COLUMBIA BIOSCIENCES**