



Item Number M1063
Item Magnesium Sulfate, Anhydrous, Powder,
CAS Number 7487-88-9
Molecular Formula MgSO₄
Molecular Weight 120.37

Assay (MgSO₄)..... Min. 99%

MAXIMUM LIMITS

Ammonium (NH₄)..... 0.005%
Arsenic (As)..... 0.0001%
Calcium (Ca)..... 0.04%
Chloride (Cl)..... 0.005%
Heavy Metals (as Pb)..... 0.001%
Iron (Fe)..... 0.001%
Loss on Ignition..... 5.0%
Manganese (Mn)..... 0.001%
Nitrate (NO₃)..... 0.005%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843