

Fighting infection – the VWR advantage

01. ERGONOMIC SEATING
SOLUTIONS

02. IDEAL FOR TECHNICAL & SCIENCE
WORK ENVIRONMENTS





Protecting front-line professionals and beyond

Now more than ever, easy and effective cleaning and disinfection of healthcare and laboratory environments, equipment and furnishings are critical. That's why Avantor offers a comprehensive collection of ergonomic seating specially constructed and designed with durable materials that stand up to more frequent and longer exposure to cleaners and disinfectants.

The VWR® lab chair line is composed of a variety of seating styles and configurations to address multidiscipline tasks and applications in technical and scientific work environments. Models include comfort-enhancing upholstered chairs and stools and polypropylene chairs highlighted by antimicrobial properties and UV inhibitors. Each design features in-demand ergonomic functionality to keep workers performing at their best and is engineered for long-term, multi-shift use.



Avantor is a globally trusted provider of mission-critical products and services enabling customers to move science forward from scientific discovery to scale-up and commercial delivery through convenient access with our premier delivery channel, VWR.

Upholstered seating: cleanability and ergonomic performance

Backed by a 10-year warranty (5 years on upholstery)



VWR® UPHOLSTERED LAB CHAIRS WITH ARMS

Description	Seat height adjustment range, cm (in.)	Weight, kg (lbs.)	Cat. No.
With casters			
Desk/low bench height	43.2 – 55.9 (17 – 22)	14.5 (31.9)	76308-748
High bench height	57.8 – 81.3 (22.75 – 32)	16.5 (36.4)	76308-740
With glides			
Desk/low bench height	43.2 – 55.9 (17 – 22)	14.2 (31.4)	76308-752
High bench height	57.8 – 81.3 (22.75 – 32)	16.3 (35.9)	76308-744



VWR® UPHOLSTERED LAB CHAIRS WITHOUT ARMS

Description	Seat height adjustment range, cm (in.)	Weight, kg (lbs.)	Cat. No.
With casters			
Desk/low bench height	43.2 – 55.9 (17 – 22)	11.7 (25.9)	76308-772
High bench height	57.8 – 81.3 (22.75 – 32)	13.8 (30.4)	76308-764
With glides			
Desk/low bench height	43.2 – 55.9 (17 – 22)	11.5 (25.4)	76308-776
High bench height	57.8 – 81.3 (22.75 – 32)	13.6 (29.9)	76308-768

Optional VWR lab chair height-adjustable arm set with hardware

Description	Adjustment range, cm (in.)	Weight, kg (lbs.)	Cat. No.
Height adjustable arms, 1 pair	6.35 (2.5)	2.9 (6.3)	76311-740

WARNING: This product can expose you to chemicals including wood dust, which is known to the state of California to cause cancer. For more information go to www.P65Warnings.ca.gov.


Seamless plastic back cover fits snugly and cleans easily, facilitating overall efficient cleaning and helping prevent contaminant buildup between back and cover.




Smooth plastic controls wipe down in an instant.



Backrest height adjustment mechanism cover wipes down easily and aids in prevention of dust, debris and contaminant infiltration of mechanism.



Chrome-finished footings with T-shaped support construction are easy to clean and disinfect.



Vinyl upholstery on seats and backs can be sanitized with EPA-certified cleaners and disinfectants. VWR 5-year upholstery warranty ensures vinyl will stand up to frequent cleaning.



High-profile glass-filled composite bases efficiently and quickly clean.



Seamless urethane armrests clean easily and prevent infiltration of contaminants.



Hooded resistance casters feature smooth tops and sides, making them easy to clean and disinfect.



Polypropylene chairs: feature tough construction and built-in protection against microbes




VWR® POLYPROPYLENE LAB CHAIRS

Description	Seat height adjustment range, cm (in.)	Weight, kg (lbs.)	Cat. No.
With casters			
Desk/low bench height	38.1 – 50.8 (15 – 20)	10.2 (22.4)	76308-760
High bench height	52.7 – 76.2 (20.75 – 30)	12.2 (26.9)	76308-756
With glides			
Desk/low bench height	38.1 – 50.8 (15 – 20)	9.9 (21.9)	76308-762
High bench height	52.7 – 76.2 (20.75 – 30)	12.0 (26.4)	76308-758


⚠ WARNING: This product can expose you to chemicals including wood dust, which is known to the state of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

Backed by a
10-year warranty

Easy-to-clean seats and backs molded of polypropylene formulation with antimicrobial properties and UV inhibitors blended in to last the life of the chair, fighting against microbes and standing up to harsh ultraviolet lighting prevalent in laboratories, helping ensure long-term performance.



Chrome-finished footrings with T-shaped support construction are easy to clean and disinfect.



Proprietary easy-to-clean stipple finish of seat and back helps prevent buildup of debris and contaminants and keeps users from sliding while seated.



High-profile glass-filled composite bases clean efficiently and quickly.



Smooth plastic controls wipe down in an instant.



Hooded resistance casters feature smooth tops and sides, making them easy to clean and disinfect.



Upholstered lab stools: easy to clean with reliable ergonomic performance



VWR® UPHOLSTERED STOOLS

Description	Seat height adjustment range, cm (in.)	Weight, kg (lbs.)	Cat. No.
With casters			
Desk/low bench height	43.2 – 55.9 (17 – 22)	6.3 (13.9)	76308-788
High bench height	57.8 – 81.3 (22.75 – 32)	8.3 (18.4)	76308-780
With glides			
Desk/low bench height	43.2 – 55.9 (17 – 22)	6.1 (13.5)	76308-792
High bench height	57.8 – 81.3 (22.75 – 32)	8.2 (18)	76308-784

⚠ WARNING: This product can expose you to chemicals including wood dust, which is known to the state of California to cause cancer. For more information go to www.P65Warnings.ca.gov.



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.