

Material No.: 0563-26  
Revision No.: 0

## Product Specification

GMP Manufactured Product

Test	Specification
Appearance (clear, colorless liquid)	Passes Test
Cyanide (CN)	Passes Test
Heavy Metals (as Pb)	$\leq 10$ ppm
Identification (by FTIR)	Passes Test
ACS – Iron (Fe)	$\leq 3$ ppm
Molarity (by titrimetry)	5.8 – 6.2
pH (1:10 dilution)	4.5 – 6.5
Ultraviolet Absorbance (1.00-cm cell vs. water) – 260 nm	$\leq 0.1$

GMP Manufactured Product  
CAUTION: For Manufacturing, processing or repackaging

Packaging Site: Paris Mfg Ctr & DC