


Anti-Adrenocorticotropin Antibody [ABT-ACTH] (A95521)


Specifications:

Name:	Anti-Adrenocorticotropin Antibody [ABT-ACTH]
Description:	Mouse monoclonal (ABT-ACTH) antibody to Adrenocorticotropin.
Applications:	IHC
Reactivity:	Human
Immunogen:	Synthetic peptide derived from human Adrenocorticotropin.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	ABT-ACTH
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Immunogen affinity purified.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline with 0.5% Bovine Serum Albumin and 0.1% Kathon CG/ICP.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:




Immunohistochemical analysis of paraffin-embedded human Pituitary using Anti-Adrenocorticotropin Antibody.




Immunohistochemical analysis of paraffin-embedded human Pituitary using Anti-Adrenocorticotropin Antibody.

Anti-Adrenocorticotropin Antibody [ABT-ACTH] (A95521)

Images continued:



Immunohistochemical analysis of paraffin-embedded human Pituitary adenoma using Anti-Adrenocorticotropin Antibody.



Immunohistochemical analysis of paraffin-embedded human Pituitary adenoma using Anti-Adrenocorticotropin Antibody.