



Item Number SO181
Item Sodium Metabisulfite, Granular, FCC
CAS Number 7681-57-4
Molecular Formula $\text{Na}_2\text{S}_2\text{O}_5$
Molecular Weight 190.11
Assay ($\text{Na}_2\text{S}_2\text{O}_5$).....90.0-100.5%

MAXIMUM LIMITS

Iron (Fe).....10 mg/kg
Lead.....2 mg/kg
Selenium (Se).....5 mg/kg

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843