



Isopropyl Alcohol AR (ACS)

Product No. 3032
Specifications current as of: Jun 8 2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Appearance (clear, colorless liquid)	Passes Test
Assay (CH ₃ CHOHCH ₃)	min. 99.5 %
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001 %
Solubility in H ₂ O	Passes Test
Titration Acid or Base (meq/g)	max. 0.0001
Water (H ₂ O)(by Karl Fischer titrn)	max. 0.2 %
Carbonyl Compounds:	
Propionaldehyde	max. 0.002 %
Acetone	max. 0.002 %

For Laboratory, Research or Manufacturing Use