

# Certificate of Quality



Part Number	MFLX95821-02
Lot Number	0211976973
Part Description	TUBING COILED FEP 1/4X3/8" 2'L
Country of Origin	USA

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: Fluorinated ethylene propylene
- The raw resin used in these products are approved under the FDA Code of Federal Regulations 21, subchapter B Part 177.1550 - for use with food and food products.
- The tubing itself has not been tested for approval by the FDA or any other testing body.
- Material is manufactured in an ISO9001-2015 Approved facility.
- All FEP tubing is manufactured with virgin resin.
- The proper cleaning of this product for food applications is the responsibility of the end-user.
- Material can be sterilized in an autoclave.
- Material is RoHs Compliant.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.