



MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System

Date of issue: 02/04/2013

Version 1.0

SECTION 1. Identification

Product identifier

Product number 105341
Product name Devarda's alloy for analysis EMSURE®

Relevant identified uses of the substance or mixture and uses advised against

Identified uses Reagent for analysis

Details of the supplier of the safety data sheet

Company EMD Millipore Corporation | 290 Concord Road, Billerica, MA 01821,
United States of America | SDS Phone Support: +1-978-715-1335 |
General Inquiries: +1-978-751-4321 | Monday to Friday, 9:00 AM to
4:00 PM Eastern Time (GMT-5)

e-mail: mm_sds@merckgroup.com

Emergency telephone 613-996-6666 CANUTEC (Canada)
+1-703-527-3887 CHEMTREC (International)
24 Hours/day; 7 Days/week

SECTION 2. Hazards identification

GHS-Labeling

Not a dangerous substance according to GHS.

Other hazards

None known.

SECTION 3. Composition/information on ingredients

CAS-No. 8049-11-4

Remarks WHMIS hazardous composition: No ingredients are hazardous
according to the CPR criteria.

SECTION 4. First aid measures

Description of first-aid measures

Inhalation

After inhalation: fresh air.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 105341
Product name Devarda's alloy for analysis EMSURE®

Version 1.0

Skin contact

After skin contact: wash off with plenty of water. Remove contaminated clothing.

Eye contact

After eye contact: rinse out with plenty of water with the eyelid held wide open. Call in ophthalmologist if necessary.

Ingestion

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

Never give anything by mouth to an unconscious person.

Most important symptoms and effects, both acute and delayed

irritant effects, Cough, Shortness of breath

Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media

For this substance/mixture no limitations of extinguishing agents are given.

Special hazards arising from the substance or mixture

Not combustible.

Ambient fire may liberate hazardous vapors.

Advice for firefighters

Special protective equipment for fire-fighters

In the event of fire, wear self-contained breathing apparatus.

SECTION 6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Avoid inhalation of vapors/aerosols or dusts. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

Environmental precautions

No special precautionary measures necessary.

Methods and materials for containment and cleaning up

Observe possible material restrictions (see sections 7 and 10).

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

SECTION 7. Handling and storage

Precautions for safe handling

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 105341
Product name Devarda's alloy for analysis EMSURE®

Version 1.0

Observe label precautions.

Conditions for safe storage, including any incompatibilities

Tightly closed. Dry.

Storage temperature: no restrictions.

SECTION 8. Exposure controls/personal protection

Exposure limit(s)

Ingredients

Basis	Value	Threshold limits	Remarks
<i>General threshold limit value for dust</i>			
CAD AB OEL	Time Weighted Average (TWA):	10 mg/m ³	Form of exposure: Total particulate.
	Time Weighted Average (TWA):	3 mg/m ³	Form of exposure: Respirable particles.
CAD BC OEL	Time Weighted Average (TWA):	10 mg/m ³	Form of exposure: Total dust.
	Time Weighted Average (TWA):	3 mg/m ³	Form of exposure: Respirable fraction.
CAD MB OEL	Time Weighted Average (TWA):	10 mg/m ³	Form of exposure: Inhalable particles.
	Time Weighted Average (TWA):	3 mg/m ³	Form of exposure: Respirable particles.
CAD ON OEL	Time Weighted Average (TWAEV):	10 mg/m ³	Form of exposure: Inhalable
	Time Weighted Average (TWAEV):	3 mg/m ³	Form of exposure: Respirable particles.
OEL (QUE)	Time Weighted Average (TWA):	10 mg/m ³	Form of exposure: Total dust.

Contains no substances with occupational exposure limit values.

Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

Hygiene measures

Change contaminated clothing. Wash hands after working with substance.

Eye/face protection

Safety glasses

Hand protection

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

Respiratory protection

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 105341
Product name Devarda's alloy for analysis EMSURE®

Version 1.0

SECTION 9. Physical and chemical properties

Physical state	solid
Color	gray
Odor	odorless
Odor Threshold	No information available.
pH	No information available.
Melting point	490 - 560 °C
Boiling point/boiling range	1663 °F (906 °C)
Flash point	No information available.
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Vapor pressure	No information available.
Relative vapor density	No information available.
Relative density	ca.5.79 g/cm ³ at 68 °F (20 °C)
Water solubility	at 68 °F (20 °C) insoluble
Partition coefficient: n-octanol/water	No information available.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	No information available.
Bulk density	ca.2,080 kg/m ³
Particle size	Particle size 0.05 - 0.20 mm

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 105341
Product name Devarda's alloy for analysis EMSURE®

Version 1.0

SECTION 10. Stability and reactivity

Reactivity

See below

Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

Possibility of hazardous reactions

no information available

Conditions to avoid

no information available

Incompatible materials

no information available

Hazardous decomposition products

no information available

SECTION 11. Toxicological information

Information on toxicological effects

Likely route of exposure

Eye contact, Skin contact, Ingestion

Specific target organ systemic toxicity - single exposure

The substance or mixture is not classified as specific target organ toxicant, single exposure.

Specific target organ systemic toxicity - repeated exposure

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Aspiration hazard

Regarding the available data the classification criteria are not fulfilled.

Carcinogenicity

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 105341
Product name Devarda's alloy for analysis EMSURE®

Version 1.0

Further information

Quantitative data on the toxicity of this product are not available.
Further toxicological data:
Hazardous properties cannot be excluded but are unlikely when the product is handled appropriately.
Possible effects:
After inhalation: Irritations of the mucous membranes, coughing, and dyspnoea.
Metal-fume fever after inhalation of large quantities.
Further data:
Handle in accordance with good industrial hygiene and safety practice.

SECTION 12. Ecological information

Ecotoxicity

No information available.

Persistence and degradability

No information available.

Bioaccumulative potential

No information available.

Mobility in soil

No information available.

Other adverse effects

SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

SECTION 14. Transport information

Land transport (DOT)

Not classified as dangerous in the meaning of transport regulations.

Air transport (IATA)

Not classified as dangerous in the meaning of transport regulations.

Sea transport (IMDG)

Not classified as dangerous in the meaning of transport regulations.

SECTION 15. Regulatory information

United States of America

Canada

WHMIS Classification

Not controlled under WHMIS (Canada).

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 105341
Product name Devarda's alloy for analysis EMSURE®

Version 1.0

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR.

Notification status

TSCA: Not Listed on TSCA inventory. For Research and Development Use only. Not For Manufacturing or Commercial Purposes.

Ingredients

Devarda alloy

DSL: This product contains one or several components that are not on the Canadian DSL nor NDSL.

Ingredients

Devarda alloy

SECTION 16. Other information

Training advice

Provide adequate information, instruction and training for operators.

Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at www.wikipedia.org.

Date of issue: 02/04/2013

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.