

VWR® GLOVE LINERS

Ideal for Critical Cleanroom
Environments

Comfortable Barrier
Wicks Moisture During
Extended Use

Low-Lint, 100%
Nylon is Durable
& Launderable





**For use under Latex,
Nitrile, PVC, and all
Disposable Gloves**

- Reduces perspiration and wicks moisture during extended use
- Prevents skin irritation without diminishing tactile sensitivity
- Limits wear and tear to extend life of gloves and conserve resources
- Low-lint, 100% nylon composition, 15 denier
- Sturdy, launderable, and reusable
- Independent users have reported 10-15 washes
- Heat-set to limit shrinkage
- Compatible with stringent cleanroom requirements
- 100% quality control inspection
- 20 pair per bag, 15 bags per case



Specifications	
Basic Weight:	21.5g/pair ± 10%
Thickness:	2mm ± 10%
Surface Resistivity:	1x10 ⁹ ± 10% ohms/square

Size	Cuff Color	Cat. No.	Pack of	Case of
Full Finger Liners				
Small	Green	94000-608	20 pr	300 pr
Medium	Orange	94000-610	20 pr	300 pr
Large	Blue	94000-612	20 pr	300 pr
X-Large	Red	94000-614	20 pr	300 pr
Half Finger Liners				
Small	Green	94000-616	20 pr	300 pr
Medium	Orange	94000-618	20 pr	300 pr
Large	Blue	94000-620	20 pr	300 pr
X-Large	Red	94000-622	20 pr	300 pr



1.800.932.5000 | vwr.com

11433-LM-LH

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2018 VWR International, LLC. All rights reserved.