



# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Revision Date 22.01.2014

Version 6.2

## SECTION 1. Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

Catalogue No.	106265
Product name	Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262
REACH Registration Number	01-2119485123-42-XXXX
CAS-No.	6131-90-4

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Pharmaceutical production, Cosmetic raw material For additional information on uses please refer to the Merck Chemicals portal ( <a href="http://www.merck-chemicals.com">www.merck-chemicals.com</a> ).
-----------------	---

### 1.3 Details of the supplier of the safety data sheet

Company	Merck KGaA * 64271 Darmstadt * Germany * Phone:+49 6151 72-0
Responsible Department	EQ-RS * e-mail: <a href="mailto:prodsafe@merckgroup.com">prodsafe@merckgroup.com</a>

### 1.4 Emergency telephone number

Please contact the regional company representation in your country.

## SECTION 2. Hazards identification

### 2.1 Classification of the substance or mixture

This substance is not classified as dangerous according to European Union legislation.

### 2.2 Label elements

#### Labelling (REGULATION (EC) No 1272/2008)

Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.

### 2.3 Other hazards

None known.

## SECTION 3. Composition/information on ingredients

### 3.1 Substance

Formula	$\text{CH}_3\text{COONa} \cdot 3 \text{H}_2\text{O}$	$\text{C}_2\text{H}_3\text{NaO}_2 \cdot 3 \text{H}_2\text{O}$ (Hill)
EC-No.	204-823-8	
Molar mass	136,08 g/mol	

For the full text of the H-Statements mentioned in this Section, see Section 16.

Remarks	No disclosure requirement according to Regulation (EC) No. 1907/2006.
---------	---

### 3.2 Mixture

not applicable

SAFETY DATA SHEET  
according to Regulation (EC) No. 1907/2006

Catalogue No.	106265
Product name	Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262

---

---

## SECTION 4. First aid measures

### 4.1 Description of first aid measures

After inhalation: fresh air.

After skin contact: wash off with plenty of water. Remove contaminated clothing.

After eye contact: rinse out with plenty of water with the eyelid held wide open. Call in ophthalmologist if necessary.

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

### 4.2 Most important symptoms and effects, both acute and delayed

irritant effects, Stomach/intestinal disorders

### 4.3 Indication of any immediate medical attention and special treatment needed

No information available.

---

## SECTION 5. Firefighting measures

### 5.1 Extinguishing media

*Suitable extinguishing media*

Water, Carbon dioxide (CO<sub>2</sub>), Foam, Dry powder

*Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### 5.2 Special hazards arising from the substance or mixture

Combustible.

Forms explosive mixtures with air on intense heating.

Development of hazardous combustion gases or vapours possible in the event of fire.

### 5.3 Advice for firefighters

*Special protective equipment for firefighters*

In the event of fire, wear self-contained breathing apparatus.

*Further information*

Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## SECTION 6. Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Avoid inhalation of dusts. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### 6.2 Environmental precautions

Do not empty into drains.

### 6.3 Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

### 6.4 Reference to other sections

Indications about waste treatment see section 13.

SAFETY DATA SHEET  
according to Regulation (EC) No. 1907/2006

Catalogue No. 106265  
Product name Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262

---

---

## SECTION 7. Handling and storage

### 7.1 Precautions for safe handling

#### *Advice on safe handling*

Observe label precautions.

#### *Hygiene measures*

Change contaminated clothing. Wash hands after working with substance.

### 7.2 Conditions for safe storage, including any incompatibilities

#### *Storage conditions*

Tightly closed. Dry.

Recommended storage temperature see product label.

### 7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

---

## SECTION 8. Exposure controls/personal protection

### 8.1 Control parameters

### 8.2 Exposure controls

#### **Engineering measures**

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

See section 7.1.

#### **Individual protection measures**

Protective clothing needs to be selected specifically for the workplace, depending on concentrations and quantities of the hazardous substances handled. The chemical resistance of the protective equipment should be enquired at the respective supplier.

#### *Eye/face protection*

Safety glasses

#### *Hand protection*

full contact:

Glove material:	Nitrile rubber
Glove thickness:	0,11 mm
Break through time:	> 480 min

splash contact:

Glove material:	Nitrile rubber
Glove thickness:	0,11 mm
Break through time:	> 480 min

The protective gloves to be used must comply with the specifications of EC Directive 89/686/EEC and the related standard EN374, for example KCL 741 Dermatril® L (full contact), KCL 741 Dermatril® L (splash contact).

The breakthrough times stated above were determined by KCL in laboratory tests acc. to EN374 with samples of the recommended glove types.

SAFETY DATA SHEET  
according to Regulation (EC) No. 1907/2006

Catalogue No.	106265
Product name	Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262

This recommendation applies only to the product stated in the safety data sheet, supplied by us and for the designated use. When dissolving in or mixing with other substances and under conditions deviating from those stated in EN374 please contact the supplier of CE-approved gloves (e.g. KCL GmbH, D-36124 Eichenzell, Internet: [www.kcl.de](http://www.kcl.de)).

*Respiratory protection*

required when dusts are generated.

Empfohlener Filtertyp: Filter P 1 (acc. to DIN 3181) for solid particles of inert substances

The entrepreneur has to ensure that maintenance, cleaning and testing of respiratory protective devices are carried out according to the instructions of the producer. These measures have to be properly documented.

**Environmental exposure controls**

Do not empty into drains.

---

**SECTION 9. Physical and chemical properties**

**9.1 Information on basic physical and chemical properties**

Form	solid
Colour	colourless
Odour	weakly of acetic acid
Odour Threshold	No information available.
pH	7,5 - 9,2 at 50 g/l 20 °C
Melting point	58 °C
Boiling point	No information available.
Flash point	> 250 °C Method: c.c. (anhydrous substance)
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Vapour pressure	at 20 °C not applicable
Relative vapour density	No information available.
Relative density	1,42 g/cm <sup>3</sup> at 20 °C

**SAFETY DATA SHEET**  
according to Regulation (EC) No. 1907/2006

Catalogue No.	106265
Product name	Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262

---

Water solubility	613 g/l at 20 °C
Partition coefficient: n-octanol/water	log Pow: -4,22 (anhydrous substance) Bioaccumulation is not expected.
Auto-ignition temperature	No information available.
Decomposition temperature	> 120 °C Elimination of water of crystallisation
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Oxidizing properties	none

**9.2 Other data**

Ignition temperature	607 °C
Bulk density	ca.900 kg/m <sup>3</sup>

---

**SECTION 10. Stability and reactivity**

**10.1 Reactivity**

Forms explosive mixtures with air on intense heating.  
A range from approx. 15 Kelvin below the flash point is to be rated as critical.  
The following applies in general to flammable organic substances and mixtures: in correspondingly fine distribution, when whirled up a dust explosion potential may generally be assumed.

**10.2 Chemical stability**

The product is chemically stable under standard ambient conditions (room temperature) .

**10.3 Possibility of hazardous reactions**

Risk of explosion with:  
nitrates  
Exothermic reaction with:  
Fluorine

**10.4 Conditions to avoid**

Strong heating.

**10.5 Incompatible materials**

no information available

**10.6 Hazardous decomposition products**

no information available

SAFETY DATA SHEET  
according to Regulation (EC) No. 1907/2006

Catalogue No.	106265
Product name	Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262

---

---

## SECTION 11. Toxicological information

### 11.1 Information on toxicological effects

#### *Acute oral toxicity*

LD50 rat: 3.530 mg/kg (anhydrous substance) (RTECS)

Symptoms: Stomach/intestinal disorders

#### *Acute inhalation toxicity*

LC50 rat: > 30 mg/l; 1 h (anhydrous substance) (RTECS)

Symptoms: Possible damages: slight mucosal irritations

#### *Acute dermal toxicity*

LD50 rabbit: > 10.000 mg/kg (anhydrous substance) (RTECS)

#### *Skin irritation*

Possible damages: slight irritation

#### *Eye irritation*

Possible damages: slight irritation

#### *Sensitisation*

This information is not available.

#### *Germ cell mutagenicity*

This information is not available.

#### *Carcinogenicity*

This information is not available.

#### *Reproductive toxicity*

This information is not available.

#### *Teratogenicity*

This information is not available.

#### *Specific target organ toxicity - single exposure*

This information is not available.

#### *Specific target organ toxicity - repeated exposure*

This information is not available.

#### *Aspiration hazard*

This information is not available.

### 11.2 Further information

Handle in accordance with good industrial hygiene and safety practice.

---

## SECTION 12. Ecological information

### 12.1 Toxicity

#### *Toxicity to fish*

LC50 *Lepomis macrochirus* (Bluegill sunfish): 5.000 mg/l; 24 h (anhydrous substance) (Lit.)

#### *Toxicity to daphnia and other aquatic invertebrates*

EC50 *Daphnia magna* (Water flea): > 1.000 mg/l; 48 h (anhydrous substance) (IUCLID)

#### *Toxicity to bacteria*

EC50 *Pseudomonas putida*: 7.200 mg/l; 18 h (anhydrous substance) (IUCLID)

### 12.2 Persistence and degradability

SAFETY DATA SHEET  
according to Regulation (EC) No. 1907/2006

Catalogue No. 106265  
Product name Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262

---

*Biodegradability*

99 %; 28 d  
OECD Test Guideline 301A  
(anhydrous substance)  
Readily biodegradable.

**12.3 Bioaccumulative potential**

*Partition coefficient: n-octanol/water*  
log Pow: -4,22  
(anhydrous substance) Bioaccumulation is not expected.

**12.4 Mobility in soil**

No information available.

**12.5 Results of PBT and vPvB assessment**

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted.

**12.6 Other adverse effects**

Discharge into the environment must be avoided.

---

**SECTION 13. Disposal considerations**

*Waste treatment methods*

See [www.retrologistik.com](http://www.retrologistik.com) for processes regarding the return of chemicals and containers, or contact us there if you have further questions.

---

**SECTION 14. Transport information**

**Land transport (ADR/RID)**

**14.1 - 14.6** Not classified as dangerous in the meaning of transport regulations.

**Inland waterway transport (ADN)**

Not relevant

**Air transport (IATA)**

**14.1 - 14.6** Not classified as dangerous in the meaning of transport regulations.

**Sea transport (IMDG)**

**14.1 - 14.6** Not classified as dangerous in the meaning of transport regulations.

**14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

Not relevant

---

**SECTION 15. Regulatory information**

**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

*EU regulations*

Major Accident Hazard 96/82/EC  
Legislation Directive 96/82/EC does not apply

Regulation (EC) No 1005/2009 on substances that deplete the ozone layer not regulated

**SAFETY DATA SHEET**  
according to Regulation (EC) No. 1907/2006

Catalogue No.	106265
Product name	Sodium acetate trihydrate extra pure Ph Eur,BP,JP,USP,FCC,E 262

---

Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC	not regulated
---	---------------

Regulation (EC) No 689/2008 concerning the export and import of dangerous chemicals	not regulated
---	---------------

Substances of very high concern (SVHC)	This product does not contain substances of very high concern above the respective regulatory limit (> 0.1 % (w/w) Regulation (EC) No 1907/2006 (REACH), Article 57).
--	---

*National legislation*

Storage class	10 - 13
---------------	---------

**15.2 Chemical Safety Assessment**

For this product a chemical safety assessment was not carried out.

---

**SECTION 16. Other information**

**Training advice**

Provide adequate information, instruction and training for operators.

**Labelling (67/548/EEC or 1999/45/EC)**

The product does not need to be labelled in accordance with EC directives or respective national laws.

**Key or legend to abbreviations and acronyms used in the safety data sheet**

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

**Regional representation**

This information is given on the authorised Safety Data Sheet for your country.

---

*The information contained herein is based on the present state of our knowledge. It characterises the product with regard to the appropriate safety precautions. It does not represent a guarantee of any properties of the product.*