



Item Number CH135
Item Cholesterol, NF
CAS Number 57-88-5
Molecular Formula $C_{27}H_{46}O$
Molecular Weight 386.65

Solubility in Alcohol..... To pass test
Melting Range..... 147°-150°C
Specific Rotation..... -34° to -38°

MAXIMUM LIMITS

Acidity..... To pass test
Loss on Drying..... 0.3%
Residue on Ignition..... 0.1%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843