



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1123-4LG, BDH1123-4LP, BDH1123-19L, BDH1123-204L, CABDH1123-4LG, CABDH1123-4LP, CABDH1123-19L, CABDH1123-204L

Catalog Description: ETHYL ACETATE

Grade: ACS

Date Current Specification Set: February 13 2007

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance	Clear, colorless liquid
Water (K.F.):	Less than 0.20%
Residue:	Less than 0.003% (0.0012 g/45mL)
GC Purity:	Greater than 99.5%
Color (APHA):	Less than 10
Titration acid:	Less than 0.9 µeq/g (Not more than 0.80 mL of 0.01N sodium hydroxide)
Substances darkened by sulfuric acid:	Passes test (No dark color at interface)

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor