

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Class 100 sterile nitrile gloves, White colour
Country of Origin	Manufactured in Malaysia
Date of Issue (yyyy-mm-dd)	2021-09-22

VWR Catalogue Number	Glove Size
76423-444	X-Small
76423-446	Small
76423-448	Medium
76423-450	Large
76423-452	X-Large
76423-454	2X-Large

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

- Water Extractable Ions Test **Pass**
- Liquid Particle Count Test **Pass**
- Total Non-Volatile Residue Test **Pass**
- Fourier Transform Infra-Red Test **Pass**
- Electrical Properties Test **Pass**
- Watertight test **Inspection level =G1; AQL 1.5**
- Visual inspection **Inspection level =G1; AQL 1.5 (critical defects)
AQL 2.5 (major defects)
AQL 4.0 (minor defects)**
- Dimensional **Inspection level = S2; AQL 4.0**
- Physical Property **Inspection level = S2; AQL 4.0**
- Glove count **Inspection level = 5 bags; Tolerance +/- 1 pair**

Product Specifications

Material	Acrylonitrile-Butadiene.
Sterilization	Product is Gamma irradiated and dose released upon ISO11137-1 recommended practices in effect at the time of validation. Products meets a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 25-40 kGy.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Recommended storage temperature is 10 °C to 40 °C. Avoid exposure to naked flames, direct sunlight or ozone source and store in a cool dry place.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only. This document has been produced electronically and is valid without a signature.