

Product Name:	VWR-HC-RFC9G
Description:	<p>VWR Catalog No. 76580-182  VWR Plus Series Refrigerator &amp; Freezer Combination  9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer)  Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 1°C – 10°C;  Freezer: -15°C – -25°C  Audible and visual high and low temperature alarms, refrigerator and freezer  2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall  Back wall evaporator cover, Refrigerator  Forced draft circulation Refrigerator  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handle, freezer; heavy duty door handle, refrigerator  Leveling legs  Stacking kit  LED interior lights  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Quality certified with 100% unit temperature testing  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 ¾ "W x 25 ½ "D x 66 ¾ "H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator and 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 269 lbs.  *Please note this unit requires 2 outlets.</p>