

# VWR®

**vwr**™  
part of avantor

# electrophoresis power supplies

01. IDEAL FOR NUCLEIC  
ACID AND PROTEIN  
APPLICATIONS

02. RUN MULTIPLE  
GELS AT ONCE

03. COMPACT SPACE  
SAVING DESIGN





### VWR® 300V AND 500V POWER SUPPLY

- Color touch screen control panel
- Compact footprint, takes up minimal lab space
- Four sets of output terminals compatible with all popular gel tanks
- Ideal for horizontal DNA gels, vertical protein gels

#### Specifications

	300V (1 volt increments)	500V (1 volt increments)
<b>Voltage Output</b>	300V (1 volt increments)	500V (1 volt increments)
<b>Current Output</b>	2000mA (5 mA increments)	400mA (1 mA increments)
<b>Output Control</b>	Constant voltage or constant current	Constant voltage or constant current
<b>Maximum Power</b>	300W	120W
<b>Output Terminals</b>	4 pairs	4 pairs
<b>Electrical Input</b>	115V / 60Hz or 230V / 50Hz	115V / 60Hz or 230V / 50Hz
<b>Weight/Dimensions</b>	1.2 kg (2.7 lbs)/ 8.9 x 12.7 x 21.6 cm (3½ x 5 x 8½")	1.14 kg (2.5 lbs) / 8.9 x 12.7 x 21.6 cm (3½ x 5 x 8½")

Description	Cat. No.
300V, 115 VAC/60Hz Input	76196-458
300V, 230 VAC/50Hz Input	76196-460
500V, 115 VAC/60Hz Input	76196-454
500V, 230 VAC/50Hz Input	76196-456

### VWR® MINI POWER SUPPLY

- Compact design, fits almost anywhere
- Simple operation — connect gel box and select voltage to start
- Selectable, 100V or 200V output
- Ideal for basic DNA separations in agarose gels
- Compatible with standard horizontal gel boxes
- Two pairs of output jacks for running two gels at one time

#### Specifications

<b>Voltage Output</b>	100V or 200V, constant
<b>Current Output</b>	Maximum 200mA
<b>Maximum Power</b>	40 Watts
<b>Output Terminals</b>	2 pairs in parallel
<b>Electrical Input</b>	115V / 60 Hz or 230V / 50 Hz
<b>Weight/Dimensions</b>	0.4 kg (0.9 lbs.)/ 16 x 11.5 x 5.7 cm (5¼ x 4½ x 2¼")
<b>Warranty</b>	2 years

Description	Cat. No.
115 VAC/60 Hz Input	76196-450
230 VAC/50 Hz Input	76196-452



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.