


## Cole-Parmer® Water Bath Comparison Guide



Specifications	WB-200 Series			WB-300 Series			WB-400 Series				
Item Numbers	78905-20 78905-21	78905-22 78905-23	78905-24 78905-25	12122-91	12122-71	12122-81	12105-84 12105-85	12105-86 12105-87	12105-88 12105-89	12105-90 12105-91	12105-92 12105-93
Capacity	6.6 L	11 L	27 L	6 L	15 L	24 L	2 L	5 L	10 L	20 L	28 L
Temperature Range	Ambient +5 to 100 °C			25 to 99.9 °C			Ambient +5 to 100 °C				
Temperature Stability	±0.2 °C			±0.5 °C			±0.2 °C				
Heater Wattage	600 W	800 W	1000 W	350 W	750 W	1400 W	120 W	360 W	1000 W	1400 W	1400 W
Low Water Sensor	—	—	—	●	●	●	—	—	—	—	—
Drain Valve	●	●	●	●	●	●	●	●	●	●	●
Programmable Timer	●	●	●	—	—	—	●	●	●	●	●
Programmable Presets	—	—	—	—	—	—	●	●	●	●	●
Overtemperature Protection	●	●	●	●	●	●	●	●	●	●	●
Delayed Start	●	●	●	—	—	—	—	—	—	—	—
Visual Alarms	●	●	●	—	—	—	●	●	●	●	●
Multiple Language Interface	—	—	—	—	—	—	●	●	●	●	●
Lid Included	●	●	●	●	●	●	●	●	●	●	●
Electrical	120 VAC, 60 Hz or 220 VAC, 50/60 Hz			230 VAC, 50 Hz			120 VAC, 60 Hz or 240 VAC, 50 Hz				
Warranty	2 years			3 years			2 years				