



## VWR® TEST TUBE RACKS

As per the manufacturer, the below product meets the following criteria:

|                                  |                                |
|----------------------------------|--------------------------------|
| <b>VWR Catalog Number</b>        | 89215-728                      |
| <b>VWR Europe Catalog Number</b> | 212-0388                       |
| <b>Lot Number</b>                | V39-040320                     |
| <b>Description</b>               | VWR Rack 16mm-36-PI Nylon Red. |
| <b>Date of Manufacture</b>       | March 04, 2020                 |
| <b>Country of Origin</b>         | India.                         |

### Quality System Compliance

Products are manufactured, inspected and found to be in compliance with product and quality specification requirements as documented in the ISO 9001:2015 & ISO13485:2016 Quality Management System. Inspection records are reviewed and signed for product release.

### QC Testing

Representative production samples are collected and inspected in accordance with current applicable standard protocols.

| Test Name   | Test Reference | Accuracy                         | Result           |
|-------------|----------------|----------------------------------|------------------|
| Autoclaving | QWI-QC-07      | General Insp. Level- I, AQL- 1.5 | PASS             |
| Dimensional | QWI-QC-23      | General Insp. Level- I, AQL- 1.5 | Within tolerance |

### Product Specifications

- Material:** 30% Glass-filled Nylon with color masterbatch.
- BSE/TSE Statement:** The above-mentioned product is manufactured without the use of Raw materials of animal origin. Therefore, BSE/TSE issue does not concern this product to the best of the manufacturer's knowledge.
- Non-Pyrogenic:** Not Applicable.
- DNase & RNase Free:** Not Applicable.
- Cytotoxicity:** Not Applicable.
- Latex-Free Statement:** No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Signed:



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** October 13, 2020