



# TRIO.BAS<sup>TM</sup>

BIOLOGICAL AIR SAMPLER


Now With IQ, OQ, PQ Documents

Microbial  
High Pressure Gas  
Sampler Kits

 **VWR**<sup>TM</sup>  
part of avantor | order on  
VWR.COM


## TRIO GAS

Designed to test the presence of microorganisms in compressed gas samples supplied from tanks and pipes under pressure.



### Each kit comes with the following:

TRIO.GAS system (stainless steel electrovalve, gas connection, stainless steel fixing system for air sampler), TRIO BAS unit (MINI or MONO), battery charger and cable, ROBUSTUS carrying case (for TRIO GAS system), aspirating head, cover head, LIGHT carrying case (for TRIO BAS unit), and fillable IQ, OQ, PQ documents



### Important Note:

If **automatically** timed samples and electronic data records are required, the MINI or MONO instrument used concurrently with the GAS system is necessary.


If **manually** timed samples and hand recorded sampling data is acceptable, the MINI or MONO instrument is not necessary. Purchase the GAS system 76196-434 or 76196-436, which includes aspirating gas chambers. Choose kit based on plate size preference. (see back page)

VWR Cat. no.

TRIO GAS System + MONO air sampler, 100 liters/min., contact plate	76196-430
TRIO GAS System + MONO air sampler, 100 liters/min., Petri plate	76196-432
TRIO GAS System + MINI air sampler, 100 liters/min., contact plate	76196-438
TRIO GAS System + MINI air sampler, 100 liters/min., Petri plate	76196-440


According to ISO 8573-7

# TRIO GAS and ASPIRATING CHAMBER



**Each kit comes with the following:**  
TRIO GAS system (stainless steel electrovalve, gas connection, stainless steel fixing system for air sampler), aspirating gas chamber, aspirating head, cover head, and ROBUSTUS carrying case (for TRIO GAS system).  
\*Timer not included.

76196-436



VWR Cat. no.

TRIO GAS System + Aspirating Chamber, contact plate

**76196-434**

TRIO GAS System + Aspirating Chamber, Petri plate

**76196-436**



800 932 5000 | [VWR.COM](http://VWR.COM)

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.

Contact your Sales Representative  
or visit [www.vwr.com](http://www.vwr.com)