



Item Number E1001
Item Disodium EDTA, Dihydrate, FCC
CAS Number 6381-92-6
Molecular Formula $C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$
Molecular Weight 372.24

Assay ($C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$).....99.0-101.0%
pH of a 1 in 100 Solution..... 4.3-4.7

MAXIMUM LIMITS

Calcium..... To pass test
Lead..... 10 mg/kg
Nitrilotriacetic Acid.....0.1%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843