



# Normal Serums

ImmunoReagents' normal serum is used as a blocking agent or as a control for immunoassays. All serum is obtained from healthy, disease free animals. Serum is lipid extracted and buffer exchanged against phosphate buffer. Our serum is supplied in lyophilized form without preservative. It is recommended to reconstitute the normal serum with DI water and add a preservative of choice, or to aliquot reconstituted serum and store frozen to prevent contamination.

Product	Part Number	VWR Cat. No.	Amount
Normal Bovine Serum	10152-192	10152-190	2 ml 10 ml
Normal Cat Serum	10152-208	10152-206	2 ml 10 ml
Normal Chicken Serum	10152-198	10152-200	2 ml 10 ml
Normal Dog Serum	10152-196	10152-194	2 ml 10 ml
Normal Donkey Serum	10152-204	75806-232	2 ml 10 ml
Normal Goat Serum	10152-212	75806-252	2 ml 10 ml
Normal Guinea Pig Serum	10152-214	10152-216	2 ml 5 ml
Normal Syrian Hamster Serum	10152-250	10152-252	2 ml 5 ml
Normal Horse Serum	10152-220	10152-218	2 ml 10 ml
Normal Human Serum	10152-224	10152-222	2 ml 10 ml
Normal Mouse Serum	10152-230	75806-254	2 ml 5 ml
Normal Pig (Swine) Serum	10152-236	10152-234	2 ml 10 ml
Normal Rabbit Serum	10152-240	10152-238	2 ml 10 ml
Normal Rat Serum	10152-242	10152-244	2 ml 5 ml

PRIMARY ANTIBODIES

SECONDARY ANTIBODIES

ANTIBODY PURIFICATION REAGENTS



BULK IMMUNOASSAY REAGENTS

ELISA & WESTERN BLOT REAGENTS

PURIFIED PROTEINS & IMMUNOGLOBULINS



1.800.932.5000  
vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

0717 Lit. No. 090004W