



Product No. 3462
Specifications current as of: Jun
8 2011

TEST	SPECIFICATION
GMP Manufactured Product	
Meets U.S.P. & F.C.C. Requirements	
Assay (CH ₃ COONa) (dried basis)	99.0 - 101.0 %
Alkalinity	max. 0.05 %
Calcium and Magnesium	Passes Test
Chloride (Cl)	max. 0.035 %
Heavy Metals (as Pb)	max. 10 ppm
Identification A	Passes Test
Identification B	Passes Test
Insoluble Matter	max. 0.05 %
Lead (Pb)	max 2 mg/kg
Loss on Drying	38.0 - 41.0 %
pH...	7.5 - 9.2
Potassium Compounds	Passes Test
Sulfate (SO ₄)	max. 0.005 %

No Class 1, 2, 3 or other solvent is likely to be present in the final product.

Bulk Food Chemical

Bulk Pharmaceutical Chemical

CAUTION: For manufacturing, processing or repacking.

For a Certificate of Lot Analysis, visit www.avantormaterials.com.