

DURAN

STAINLESS STEEL BOTTLES

Ideal for storage and shipping applications where glass is not applicable, due to the risk of breakage or chemical incompatibility. Manufactured from corrosion resistant AISI Type 316 L (1.4404) stainless steel, and hygienically constructed with all welded seams. The bottle has a smooth inner surface finish that corresponds to IIIc (DIN 17441).

Lab and Storage Bottle:

- Supplied without a cap (compatible with all the GL45 caps)
- Hazardous materials laboratory container for storage of liquids, intermediates, or solid product
- Storage of precious materials, such as high purity fine chemicals, pharmaceutical or cosmetic products
- Storage of light sensitive materials

UN-Certified Bottle adds:

- Supplied with GL45 stainless steel cap and PTFE sealing disc
- Certified to UN standards for the carriage of HAZMAT Class 3: Flammable Liquids; classified as dangerous goods in UN packing Group I: (great danger) Groups II (medium danger) and III (low danger)

Stainless Steel Caps

- GL45 thread



Also Available...

Contact us for other bottles and accessories.



DURAN

LABORATORY GLASS BOTTLES

Proven through decades of practical use, DURAN® products set benchmarks where resistance and safety are concerned. Particularly when compared with plastic, the individual material properties of DURAN® offer clear advantages. Choose original DURAN® laboratory glass bottles and accessories because top quality tolerates no compromises.

GL Thread Bottles:

- Plain, safety coated, amber, or pressure versions available
- 10mL through 20L
- Graduated
- Large labeling field for easy marking
- GL25, GL32 & GL45 threads available
- Available with or without cap and pouring ring

Wide Neck Bottles:

- Plain, safety coated, or amber versions available
- Larger GL80 thread opening
- 250mL through 20L

7622 Caps and Pouring Rings

- GL threads
- Available in polybutylene, polypropylene, TpCh260 (similar to PFA)
- Various colors for safe identification
- Available PTFE lined



Also Available...

Contact us for other bottles, caps, and accessories.



avantor[™]
delivered by **VWR**[™]

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.