



Material	92368-700
Material description	Toluene
Grade	ACS
CAS Number	108-88-3
Molecular formula	C ₇ H ₈
Molecular mass	92,14

Characteristics

Specifications

Assay (GC, corrected for water)	Min. 99,5 %
Substances darkened by sulphuric acid	Passes test
Color (Pt-Co)	Max. 10 APHA
Residue after evaporation	Max. 0,001 %
Sulphur compounds (as S)	Max. 0,003 %
Water	Max. 0,03 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada