



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1111-4LG , BDH1111-19L, BDH1111-204L, BDH1111-4LG, CABDH1111-19L, CABDH1111-204L

Catalog Description: Cyclohexane

Grade: ACS

Date Current Specification Set: August 15 2007

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance:	Clear
Water (K.F.):	Less than 0.02%
Residue:	Less than 0.002% (0.0020 g/129 mL)
GC Purity:	Greater than 99.0%
Color (APHA):	Less than 10
Substances darkened by sulfuric acid:	Passes test (Neither the cyclohexane nor acid layer should be darkened.)

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor