

T 50 Digital ULTRA-TURRAX® IKA® Disperser

Powerful dispersing instrument for volumes from 0.25 - 30 l (H₂O) with digital speed display. It offers a wide speed range from 500 – 10,000 rpm that enables users to work at high circumferential speeds even with small rotor diameters. A broad choice of dispersing elements guarantees a wide range of applications. Reproducible operations due to constant speed even with changes in viscosities.

- Digital speed display
- Electronic speed control
- Electronic overload protection
- Stainless steel dispersing elements can be cleaned quickly and easily
- Smooth start
- Error code display



Motor protection against overload



Digital display for precise monitoring of set and actual speeds



Wide selection of dispersing tools to suit your application



Rotating knob for adjusting the speed

Ordering information: Three-year manufacturer's warranty. Unit includes storage case, power cord and tool kit. Dispersing elements, mounting clamps, stands and other accessories are sold separately.