

Certificate of Quality



Part Number	MFLX30547-02
Lot Number	A74762
Part Description	SANITARY GASKETS 3/4IN TRI-CLAMP 10/PACK
Expiration Date	01JAN2026
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- That all the materials of construction are USP Class VI, 21 CFR and European Communities EMA/410/01 Rev 3 compliant.
- The manufacturer has indicated that to the best of their knowledge this product does not contain animal fats or by-products including Bovine, Ovine or Caprine Tissues.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.
- This product is not approved or intended for, and should not be used for medical, clinical, surgical or other patient oriented applications.
- *Note* - we recommend visual inspection of gaskets for discoloration, hardening and deformation after three years and before use.
- Products stored in their original packaging in a dry, cool environment away from direct sun light and artificial light should remain in optimal condition for 10 years.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.