

Certificate of Conformity

Product Information

Product Number: 28903471

Product Name: illustra GFX PCR DNA and Gel Band Purification Kit

Pack Size: 250 purifications

Lot Number: 16959251

Expiry Date: 25 August 2020

Contents

28929636 Capture buffer type 3, Lot : 16959251

1 x 140 ml

9602CC Wash buffer type 1*, Lot : 16959251

1 x 62.5 ml

380180 Elution buffer type 4, Lot : 16959251

1 x 28 ml

9602DB Elution buffer type 6, Lot : 16959251

1 x 28 ml

28929969 illustra™ GFX™ MicroSpin™ columns,

Lot 312526-10

5 x 50 pack

29042200 Collection tubes: 5 x 50 pack

* Absolute ethanol provided by the researcher is added immediately prior to use. Once ethanol is added to the Wash buffer, store UPRIGHT and AIRTIGHT

Functional assay

GFX PCR DNA and Gel Band Purification Kit is tested for the ability to quantitatively isolate a 910 base-pair PCR product from solution and from an agarose gel band.

Functional test

910 bp PCR fragment purified and analysed on an agarose gel for size and yield.

Conformance & Quality systems statement

This is to certify that this product conforms to GE Healthcare Life Sciences specifications. All products are manufactured via a Quality System certified to **BS EN ISO9001:2008** and tested in accordance with documented quality procedures and approved as a result of meeting the required specification

Electronic signature

This document has been electronically produced and is valid without a signature.

Version AL

www.gelifesciences.com