

## T 25 Digital ULTRA-TURRAX® IKA® Disperser

### red**dot** design award winner 2012

High-performance dispersing instrument for volumes from 1 – 2000 ml (H<sub>2</sub>O) with digital speed display. It offers a wide speed range from 3000 – 25,000 rpm that enables users to work at high circumferential speeds even with small rotor diameters. A broad choice of dispersing elements adds versatility. Applications ranges from homogenizing waste water samples, use in laboratory reactors, dispersion tasks under vacuum / pressure and sample preparation in medical diagnostics.

- Digital speed display
- Electronic speed control
- Electronic overload protection
- Stainless steel dispersing elements can be cleaned quickly and easily
- Plastic disposable dispersing elements are available in two sizes
- Error code display
- Quiet operation



Motor protection against overload



Digital display for precise monitoring of set and actual speeds



Wide selection of dispersing tools to suit your application



Rotating knob for adjusting the speed

Ordering information: Three-year manufacturer's warranty. Unit includes storage case, power cord and tool kit. Dispersing elements, mounting clamps, stands and other accessories are sold separately.