

Product: 9" (3mil- S Series) Powder Free Nitrile Examination Glove Finger Textured Gloves (M-3.5g)

PRODUCT INFORMATION

1. Material : Nitrile butadiene rubber (S - Series)

2. Particulate Test :

Powder Residue Glove
Max 2 mg/glove

3. Color : Baby Blue

4. Design and Feature : Ambidextrous, straight fingers, finger textured and beaded.

5. Bead Diameter : Min 0.5mm

6. Dimensions :

Glove size	Glove Palm width (mm) Specification	Glove Length (mm) Specification
XS	70 ± 10	Min 230
S	80 ± 10	
M	95 ± 10	
L	110 ± 10	
XL	120 ± 10	

7. Thickness :

Location of thickness measurements	*Single wall (mm) - Specification	*Single wall (mm) - Typical Value
Finger (at 13±3mm from fingertip)	Min 0.07	0.07 - 0.10
Palm (at center of palm)	Min 0.06	0.06-0.07
Cuff (at 25±5mm from cuff edge)	min 0.05	0.05 -0.06

Remark: Typical Value for reference only

8. Physical Properties

Tensile Strength (MPa)	Before aging	After Aging
Specification	min 14 Mpa	min 14 Mpa






Ultimate Elongation (%)	Before aging	After Aging
Specification	Min 500	Min 400

9. Glove Weight (Reference Only) :

Glove size	Weight (gm)
XS	2.9 ± 0.2
S	3.2 ± 0.2
M	3.5 ± 0.2
L	3.8 ± 0.2
XL	4.1 ± 0.2

10. Certificate of Release : A Certificate of Release is available upon request

11. Quality Assurance : Meet or exceed ASTM D6319 on water leak & dimension.
 Manufactured under GMP, ISO 9001, ISO 13485 & MDSAP

Prepare by: QRM (Noor Rozinaini)	Reviewed by: QA (Poppy Farah Rossa)	Approved by: Process & Technical (Husni)	Approved by: R & D (Prabhakaran Kumaresan)	Approved by: Chief Operating Officer (Hew Yew Fook)
				

CONTROLLED COPY