

# Certificate of Quality



Part Number	<b>MFLX31809-42</b>
Lot Number	<b>28962273</b>
Part Description	<b>ELBOW 90-DEG 1IN LADISH CLAMP</b>
Expiration Date	<b>08OCT2061</b>
Country of Origin	<b>USA</b>

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: PVDF
- Material grade: Sani-Pro® K PVDF.
- Material meets 21 CFR 177.2510 requirements.
- Material meets 3-A requirements for Sanitary Standards for Multiple-Use Plastic Materials Used as Product Contact Surfaces for Dairy Equipment.
- Material meets USDA requirements or use in processing or storage areas in contact with meat or poultry food products.
- Material meets USP Class VI requirements.
- No animal derived ingredients are used in the manufacturing of Sani-Pro® K PVDF material.
- Shelf life: 40 years.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.