

# VWR® glassware washers

- 01. HIGH WATER TEMPERATURE PROVIDES EFFECTIVE REMOVAL OF OIL RESIDUE
- 02. SPINDLE RACK READY UNITS FOR CLEANING INSIDE NARROW NECK GLASSWARE
- 03. UNDER COUNTER, MOBILE AND VERTICAL SPACE-SAVER MODELS AVAILABLE





Undercounter model



Freestanding model



Vertical space saver model

## Easily remove stubborn oils from inside glassware

### VWR GLASSWARE WASHERS

Removing sticky, waxy and oily residue from glassware is a challenge especially when the residue is inside narrow necked glassware such as round bottom flasks. VWR Glassware Washers are up for the challenge. The High 70°C (158°F) water temperature, along with an appropriate detergent, effectively dissolves stubborn oils leaving your glassware sparkling clean.

VWR Glassware Washers come with top and bottom stainless-steel open racks to accommodate beakers and other wide mouth glassware. Spindle rack ready models are available for thorough cleaning inside flasks and other glassware, and additional baskets, insert racks, and spindle racks are available to maximize space, improve efficiency, and protect specialty glassware.

#### Features include:

- 6 Cycle Selections to let you select optimal temperature and time, and include a customizable timed option
- Enhanced dual fan drying system helps expel moisture for faster drying
- Noise buffering insulation for an undisturbed work environment
- Three separate filters to protect pumps from debris
- Optional Deionized (DI) Rinse models use separate pump for facility provided DI water final rinse

Maximize results with spindle racks to effectively remove oils inside flasks and other glassware. Baskets and inserts keep glassware organized and in place.

Learn more at [vwr.com](http://vwr.com) keyword: **VWR Glassware Washer**

Description	Cat. No.
Undercounter washer	10002-850
Undercounter washer with DI rinse	10002-854
Undercounter washer with DI rinse, spindle rack ready	10002-858
Freestanding and mobile washer	10002-852
Freestanding and mobile washer with DI rinse	10002-856
Freestanding and mobile washer with DI rinse, spindle rack ready	10002-860
Vertical stacked washer with top DI rinse/standard rack and bottom spindle rack ready	10002-862
Vertical stacked washer with top and bottom DI rinse and spindle rack ready	10002-864
Vertical stacked washer, standard rack for use with tap water	10002-866



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2021 Avantor, Inc. All rights reserved.