



Item Number	SO187
Item	Monobasic Sodium Phosphate, Anhydrous,
CAS Number	7558-80-7
Molecular Formula	NaH ₂ PO ₄
Molecular Weight	119.98
	Assay (NaH ₂ PO ₄ ; Anhydrous Basis).....98.0-103.0%
	pH of a 5% Solution..... 4.1-4.5
	MAXIMUM LIMITS
	Water..... 2.0%
	Insoluble Substances..... 0.2%
	Chloride (Cl)..... 0.014%
	Sulfate (SO ₄).....0.15%
	Aluminum, Calcium and Related Elements..... To pass test
	Arsenic (As)..... 8 ppm
	Heavy Metals (as Pb)..... 0.002%
	Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters:	East Coast Facility:
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	