

Certificate of Quality



| | |
|-------------------|----------------------------------|
| Part Number | MFLX96412-16 |
| Lot Number | 122021001 |
| Part Description | TUBING MFLEX FDA VITON NO.16 25' |
| Country of Origin | USA |

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: EE97360B
- This product was manufactured in a facility registered to AS9100 with ISO9001.
- All ingredients used in these products are listed under the FDA code of Federal Regulations 21, subchapter B, Part 177.2600-Rubber Articles intended for repeated use in contact with food (food grade).
- Free of any materials of animal derivatives
- Free of DEHP and Bisphenol A
- Shelf Life is 20 years from the manufacturer date.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.