

## Product Specification

Meets U.S.P Requirements, GMP Manufactured Product

Test	Specification
USP – Assay (ZnO) (ignited basis)	98.0 – 102.0 %
USP – Identification A	Passes Test
USP – Identification B	Passes Test
USP – Alkalinity	Passes Test
USP – Loss on Ignition	≤ 1.0 %
USP – Carbonate and Color of Solution	Passes Test
USP – Arsenic (As)	≤ 6 ppm
USP – Lead (Pb)	Passes Test
USP – Iron and Other Heavy Metals	Passes Test
Solubility in Dilute Sulfuric Acid	Passes Test
Lead (Pb)	≤ 0.0010 %
Cadmium (Cd)	≤ 0.0015 %
Iron (Fe)	≤ 0.0020 %

GMP Manufactured Product

Bulk Pharmaceutical Chemical

CAUTION: For Manufacturing, processing or repackaging

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Metallic Residues: No metal catalysts or metal reagents, as defined by EMA Guideline EMEA/CHMP/SWP/4446/2000 , are used in the production of this material.

Packaging Site: Paris Mfg Ctr & DC