

Catalog Number: BDH9312

Catalog Description: BDH COPPER SULFATE CRYST

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
-----------------	----------------

Calcium	<= 0.005 %
Chloride	<= 0.001 %
Insolubles	<= 0.005 %
Iron	<= 0.003 %
Nickel	<= 0.005 %
Nitrogen Compounds (as N)	<= 0.002 %
Potassium	<= 0.01 %
Purity	98.0 - 102.0 %
Sodium	<= 0.02 %

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager