

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Untapered Connectors
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2021-10-07

VWR Catalogue Number	Fits tubing I.D.	Connector
46600-056	4 mm (⁵ / ₃₂ "	T Connector
46600-058	6 mm (⁷ / ₃₂ "	T Connector
46600-148	8 mm (⁵ / ₁₆ "	T Connector
46600-060	10 mm (³ / ₈ "	T Connector
46600-062	12 mm (¹⁵ / ₃₂ "	T Connector
46600-086	14 mm (⁹ / ₁₆ "	T Connector
46600-088	16 mm (⁵ / ₈ "	T Connector
46600-122	4 mm (⁵ / ₃₂ "	L Connector
46600-124	6 mm (⁷ / ₃₂ "	L Connector
46600-126	8 mm (⁵ / ₁₆ "	L Connector
46600-128	9.5 mm (³ / ₈ "	L Connector
46600-136	4 mm (⁵ / ₃₂ "	4-Way Connector
46600-138	6 mm (⁷ / ₃₂ "	4-Way Connector
46600-140	8 mm (⁵ / ₁₆ "	4-Way Connector
46600-144	12 mm (¹⁵ / ₃₂ "	4-Way Connector
46600-146	14 mm (⁹ / ₁₆ "	4-Way Connector
37001-048	16 mm (⁵ / ₈ "	4-Way Connector
46600-064	4 mm (⁵ / ₃₂ "	Y Connector
46600-066	6 mm (⁷ / ₃₂ "	Y Connector
46600-068	8 mm (⁵ / ₁₆ "	Y Connector
46600-070	9.5 mm (³ / ₈ "	Y Connector
46600-072	12 mm (¹⁵ / ₃₂ "	Y Connector

46600-092	16 mm (⁵ / ₈ "	Y Connector
-----------	---------------------------------------	-------------

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Polypropylene; Non-ADI compliant.
Sterilization	Not applicable.
ATP Assay	Not applicable.
Melamine Statement	This product does not contain Melamine.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di(2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store in a dry place. Do not put more than 8 boxes one on the other and do not expose to temperature above 50° or below – 30°.
Shelf Life	60 months from production date written on the label.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.