

**Product Specifications**

<b>VWR Part Number:</b>	<b>IR2423</b>
<b>ImmunoReagents, Inc. Part Number:</b>	SP-128-VX2
<b>Product:</b>	NORMAL CANINE SERUM 2ML
<b>Concentration:</b>	60.0 mg/ml (Bradford, IgG standards)
<b>Amount:</b>	2.0 ml
<b>Conjugate:</b>	Unconjugated
<b>Form:</b>	Lyophilized
<b>Purification:</b>	Lipid extracted and dialyzed against 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Purity:</b>	N/A
<b>Host:</b>	Canine
<b>Immunogen:</b>	N/A
<b>Buffer:</b>	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Storage:</b>	Store lyophilized material at 2-8 °C. For long term storage after reconstitution, prepare small aliquots and store at -20 °C. For storage at 2-8 °C, add a preservative to prevent growth of bacteria.
<b>Shelf Life:</b>	1 year from date of receipt. Prepare working dilution prior to use and then discard.
<b>Specificity:</b>	Non-immune
<b>Cross Reactivity:</b>	N/A
<b>Country of Origin:</b>	Normal Dog Serum was obtained from healthy animals of Chinese origin.
<b>Applications:</b>	This product is used as a blocking reagent or control for most immunoassay applications.
<b>Disclaimer:</b>	For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.