



Quality Certificate

VWR® SPECIMEN CONTAINERS, POLYPROPYLENE WITH POLYETHYLENE CAPS

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	25384-148
Lot Number:	121447KS
Description:	VWR CONTAINER NS 4.5OZ CS500
Manufacture Date:	01/01/2018
Country of Origin:	USA

The product has been inspected and released in accordance with Plant Standard Production Methods and Process Quality Assurance Specifications. All resins used in the processing of these products have been approved and released.

Materials of Construction: Polypropylene

Non-Pyrogenic: Product has not been treated for pyrogens.

DNase & RNase Free: Product has not been treated for DNase or RNase agents.

Cytotoxicity: N/A

Latex Free: Not made with natural rubber latex.

BSE/TSE: Resins used to mold this product do not contain products of bovine origin.

Signed :

Linda K. Shine
Manager - Supplier Qualification
Global Compliance
VWR International LLC

Date: 05/04/2018