

Certificate of Quality



Part Number	MFLX96105-09
Lot Number	H033L8L001
Part Description	TUBING PHARMA-50 3/8 X 5/8 50FT
Expiration Date	26JUL2026
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: Silastic Platinum-cured
- All ingredients used in these products are approved under the FDA code of Federal Regulations 21, Subchapter B, Part 177.2600 - Rubber Articles intended for repeated use in contact with food (Food-Grade).
- The tubing itself has not been tested for approval by the FDA or any other testing body.
- The elastomer meets USP Class VI requirements summary of health data.
- Complies with the criteria in 3-A Sanitary Standards for Multiple-use Rubber and Rubber-like Materials used as Product Contact Surfaces in Dairy Equipment, Number 18-02.
- The product meets 3-A Sanitary Plastics Standards number 20-18 -criteria for handling mild and milk products.
- This product has been manufactured in accordance with the principles of Good Manufacturing Practices.
- These products have been manufactured in an ISO-9001 registered facility.
- A lot control system is used to assure traceability to all raw materials.
- These products have been tested and certified to meet the requirements of the European Pharmacopoeia (EP) Monograph 3.1.9 Silicone Elastomer for Closures and Tubing.
- The proper cleaning of this product for food applications is the responsibility of the end-user.
- This tubing product is made from Platinum catalyzed silicone.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.