

T 18 Digital ULTRA-TURRAX® IKA® Disperser

Competitively priced dispersing instrument for volumes from 1 – 1500 ml (H₂O) with digital speed display. It offers a wide speed range from 3000 – 25,000 rpm that enables users to work at high circumferential speeds even with small rotor diameters. A broad choice of dispersing elements guarantees a wide range of applications.

- Digital speed display
- Electronic speed control
- Electronic overload protection
- Stainless steel dispersing elements can be cleaned quickly and easily
- Plastic disposable dispersing elements are available in two sizes
- Error code display
- Quiet operation



Motor protection against overload



Digital display for precise monitoring of set and actual speeds



Rotating knob for adjusting the speed

Ordering information: Three-year manufacturer's warranty. Unit includes storage case, tool kit and power cord. Dispersing elements, mounting clamps, stands and other accessories are sold separately.