



**Item Number** S1228  
**Item** Sodium Carbonate, Anhydrous, FCC  
**CAS Number** 497-19-8  
**Molecular Formula** Na<sub>2</sub>CO<sub>3</sub>  
**Molecular Weight** 105.99

Assay (Na<sub>2</sub>CO<sub>3</sub>; Dry Basis).....99.5-100.5%

MAXIMUM LIMITS

Lead..... 4 mg/kg  
Loss on Drying..... 1%

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843