

Special Vivaspin® Anniversary Offer

Turbo Charge Your Concentration

Vivaspin® Family

Vivaspin® is the optimal solution for any concentration or buffer exchange application!

The Vivaspin® Turbo optimized design, its sleek internal profile, ensures maximum process speeds right the way down to the last few micro liters. The UV joining technology allows for a smooth joint transition between membrane and plastic housing, allowing all of your valuable sample to be collected into the unique pipette friendly dead-stop pocket.

Vivaspin is available in a broad range of MWCOs: 3,000 - 100,000 kDa.

Key Features & Benefits

- High-speed concentration with high recovery rate
- Highest flow-through speed and maximum retention with a broad range of volume capacities and MWCOs
- Patented angular dead-stop technology allowing easy recovery and prevents the sample from drying out
- High chemical compatibility



A \$150 Value!

Buy 3 and Get 1 Free

Buy 3 of any Vivaspin® 20- Polyethersulfone 12 pack, Vivaspin® Turbo 4-Polyethersulfone 25 pack or Vivaspin® Turbo 15-Polyethersulfone in the table below and receive one pack of equal or lesser value free.

* Special offer valid only for Canada from January 1st, 2021 until June 30th, 2021. Cannot be combined with any other offer.



VWR Cat. no.	Sartorius Part	Description
CA76409-098	VS2091	Vivaspin 20 - 3,000 MWCO PES
CA76409-070	VS2011	Vivaspin 20 - 5,000 MWCO PES
CA76409-066	VS2001	Vivaspin 20 - 10,000 MWCO PES
CA76409-074	VS2021	Vivaspin 20 - 30,000 MWCO PES
CA76409-078	VS2031	Vivaspin 20 - 50,000 MWCO PES
CA76409-082	VS2041	Vivaspin 20 - 100,000 MWCO PES
CA76409-086	VS2051	Vivaspin 20 - 300,000 MWCO PES
CA76409-090	VS2061	Vivaspin 20 - 1,000,000 MWCO PES
CA76409-094	VS2071	Vivaspin 20 - 0,2µm PES
CA76408-986	VS04T91	Vivaspin Turbo 4 - 3 kDa, PES
CA76408-970	VS04T11	Vivaspin Turbo 4 - 5 kDa, PES
CA76408-966	VS04T01	Vivaspin Turbo 4 - 10 kDa, PES
CA76408-974	VS04T21	Vivaspin Turbo 4 - 30 kDa, PES
CA76408-978	VS04T31	Vivaspin Turbo 4 - 50 kDa, PES
CA76408-982	VS04T41	Vivaspin Turbo 4 - 100 kDa, PES
CA76409-062	VS15T91	Vivaspin Turbo 15 - 3 kDa MWCO, PES
CA76409-046	VS15T11	Vivaspin Turbo 15 - 5 kDa MWCO, PES
CA76409-042	VS15T01	Vivaspin Turbo 15 - 10kDa MWCO, PES
CA76409-050	VS15T21	Vivaspin Turbo 15 - 30 kDa, PES
CA76409-054	VS15T31	Vivaspin Turbo 15 - 50 kDa, PES
CA76409-058	VS15T41	Vivaspin Turbo 15 - 100 kDa, PES

Volume	100 - 500 µm	0.4 - 2 ml	0.4 - 2 ml	2 - 4 ml	2 - 6 ml	6 - 15 ml	2 - 15 ml	6 - 20 ml
Vivaspin 500	Vivaspin 2	Vivaspin 2	Vivaspin Turbo 4	Vivaspin 6	Vivaspin Turbo 15	Vivaspin 15R	Vivaspin 20	
Membranes	PES	PES	Hydrosart (RC)	PES	PES	PES	Hydrosart (RC)	PES
Cut off: Needs to be ~3x smaller than the protein size to concentrate	3 kDa 5 kDa 10 kDa 30 kDa 50 kDa 100 kDa 300 kDa 1,000 kDa 0.2 µm	3 kDa 5 kDa 10 kDa 30 kDa 50 kDa 100 kDa 300 kDa 1,000 kDa 0.2 µm	2 kDa 5 kDa 10 kDa 30 kDa	3 kDa 5 kDa 10 kDa 30 kDa 50 kDa 100 kDa	3 kDa 5 kDa 10 kDa 30 kDa 50 kDa 100 kDa 300 kDa 1,000 kDa 0.2 µm	3 kDa 5 kDa 10 kDa 30 kDa 50 kDa 100 kDa	2 kDa 5 kDa 10 kDa 30 kDa	3 kDa 5 kDa 10 kDa 30 kDa 50 kDa 100 kDa 300 kDa 1,000 kDa 0.2 µm

To purchase your 3 cases of Vivaspin concentrators, visit vwr.com, contact your VWR Sales Rep, or call 1.800.932.5000 between January 1st and June 30th, 2021. To redeem for your free product, visit vwr.com/promotions and search for Promo Code 5587. Fill in the online form before July 15th, 2021. Allow up to 8 weeks to receive your free product.



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in Canadian dollars unless otherwise noted. Offers valid in Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2021 Avantor, Inc. All rights reserved.

0321 Lit. No. 190131WREV Promo Code 5587